

NEW PRODUCT ANNOUNCEMENT



safematic

RC-010

Remote control software for operating safematic
CCU-010 coating instruments.



Extend the application range of Safematic CCU-010 coating instruments with the Windows-based **RC-010** remote control software. With this tool you have the possibility to remotely run, monitor and analyse your coating process. You can easily create and run your own recipes according to your needs to achieve reproducible results. The graphical user interface facilitates easy setup and selection of recipes. The RC-010 tool also comes ready with predefined recipes with optimal settings for the most common metal sputtering targets as well as carbon.

It couldn't be easier: Connect your laptop or computer to the CCU-010 coating system via a USB cable and get started.

Use RC-010 for:

- Glove box applications
- Education / training
- Remote monitoring

