

Safematic LC-006 *LARGE CHAMBER*

6-inch chamber for larger samples and wafers

Area of application:

The LC-006 is an enlarged version of the current glass chamber. This chamber is made of robust aluminum with a wall thickness of 8 mm, providing high structural stability. Thanks to this solid design, no additional implosion protection is required.

With a diameter of 6.2 inches (157.5 mm), the LC-006 allows the coating of larger samples of various types, significantly expanding the system's application range.



Features:

- Retrofittable due to modular design
- Simple assembly
- Samples up to 6.2 inches (157.5mm) in diameter
- Included wafer holder
- Compatible with all options

