


# **Safematic TS-010** *THERMAL SHIELD*

## Thermal shield for gentle sample treatment

### Area of application:

The TS-010 thermal shield is an option for the CT-010 carbon thread head. The option can be used when thermal sensitive samples need to be coated. Due to its shielding effect, the samples are not exposed to direct thermal radiation and are treated more gently. The TS-010 can also help to apply thinner layers of carbon.



### Features:

- Retrofittable due to modular design
- Simple assembly
- Less thermal stress for the samples
- Enables thinner layers

